WARNING: The *hpf*ast procedure section must be read before using a test! Before performing testing all operators MUST read and be familiar with this information.

COMPLEXITY: WAIVED

Any modification by the laboratory to the test system or FDA approved test system instructions will result in the test no longer meeting the requirements for waived category.

1. NAME: *hpf*ast

Rx Only

IVD

2. INTENDED USE


*hpf*ast detects the urease enzyme for the presumptive identification of *Helicobacter pylori* in gastric mucosal biopsies. It is intended for *in vitro* diagnostic use only. IN CANADA FOR LABORATORY USE ONLY

3. SUMMARY

Helicobacter pylori has now been shown to be the causative agent in most instances of acute type B gastritis (1). And, it is now known that, in the absence of aspirin, non-steroidal anti-inflammatory drugs or hypersecretory states, the bacteria is directly related to peptic ulcer disease such as duodenal and benign gastric ulcers (2). There is epidemiological data linking the presence of *H. pylori* with gastric cancer and low-grade mucosal associated lymphoid tissue (MALT) lymphoma (3,4). The eradication of *H. pylori* gastritis by antibiotics has been shown to cure peptic ulcers and prevent recurrence. (4,5).

4. PRINCIPLE

A gastric mucosal endoscopic biopsy is placed in an agar gel containing urea, two pH dye indicators: bromthymol blue and methyl red, and a paraben preservative. *hpf*ast does not contain any toxic azide preservatives.

The urease in *H. pylori* converts the urea to ammonia which raises the pH and changes the agar color, indicating a positive test.

5. REAGENTS

5.1 Storage and Stability

*hpf*ast should be stored at 2-8°C in the boxed container. *hpf*ast is stable and provides accurate results for up to 12 months. Do not use the test if there is not a light peach or yellow color or if the test is past the expiration date. Each *hpf*ast contains an expiration date and lot number to permit tracking. Do not freeze.

5.2 Disposal

Specimens and test devices must be disposed of according to appropriate guidelines for biohazardous waste.

6. SPECIMEN COLLECTION AND PREPARATION

6.1 Patient Preparation

Antibiotics and bismuth, when used in suboptimal therapy, can suppress but not eliminate the organism. Therefore, the patient should not have used these agents for several weeks prior to the test. Following inadequate therapy, the organism may regrow in a patchy manner and may not be detected by random biopsies. The patient should not have received proton pump inhibitors such as omeprazole or lansoprazole, as these drugs have been shown to have a bacteriocidal effect on the organism prior to the test.

6.2 Biopsy Location

Because of the patchy nature of *H. pylori* infection, it is generally recommended that two specimens be obtained from each patient, one from the immediate prepyloric area and a second from the proximal antrum on the lesser curvature (6). Biopsies should be obtained from normal appearing mucosa as eroded or ulcerated mucosa may have fewer organisms.

7. PROCEDURE

STEP 1: The tray with the *hpf*ast agar should be evaluated before the test. The agar should have a light peach or yellow color. A very light green color may indicate contamination and the subsequent result should be interpreted cautiously. A tray with a dark green or blue agar color should be discarded. Allow the tray to reach room temperature before use.

STEP 2: Free the specially designed biopsy pick from the tray by squeezing the tray from both sides. The biopsy pick is designed at one end with a needle tip and spatula to remove the specimen from the forceps.

STEP 3: Peel back the label and place the specimen on the top of the agar. Use the reverse end of the biopsy pick to insert the specimen deep into the agar gel. The specimen should be exposed to as much agar as possible to enhance the likelihood of a positive reaction.

STEP 4: Resal the test by firmly pressing the label back in place. Failure to reseat the test may result in **evaporation of contents**. The patient name, date and time should be noted on the label. A warming plate or small incubator heater (35-40°C) may increase the rapidity of the reaction.

8. INTERPRETATION OF RESULTS

8.1 Suggested reading times for *hpf*ast are 15, 30 and 60 minutes, 4 hours, and 24 hours after biopsy insertion. Once a positive reaction has been obtained, no further readings are required.

8.2 The dyes, methyl red and bromthymol blue, have been selected to turn various shades of green and then blue as the pH of the medium rises.

| | | |
|-----------|--------------------|-----------------------|
| pH | Color | Interpretation |
| 5.2 | yellow/light peach | negative |
| 5.8 | greenish-yellow | negative |
| 6.5 | dark green | positive (slightly) |
| 7.2 | aqua blue | positive (moderately) |
| 7.5 | dark blue | positive (markedly) |

8.3 **A yellow/light peach or greenish-yellow color** indicates the absence of the urease enzyme and is a **negative test**. The final reading on a negative test should be taken at not more than 24 hours.

8.4 **A dark green reaction within 4 hours is a positive** (presumptive diagnosis) and will turn dark blue within 24 hours. Dark green reactions at 24 hours are rare and should be interpreted in conjunction with histological or other suitable criteria.

8.5 **A blue color indicates the presence of the urease enzyme and is positive** for *H. pylori* (presumptive diagnosis). However, for the definitive diagnosis of *H. pylori*, histological examination with a suitable stain, such as Giemsa, Warthin-Starry, or Genta stain is required.

8.6 Quality Control Procedure

8.6.1 POSITIVE CONTROL: A positive control test preparation is provided with each shipment of *hpf*ast. The positive control unit dispenses approximately 30-40 microliters per drop of the urease-glycerol preparation (Sigma Catalog # U1875).

8.6.2 POSITIVE CONTROL TESTING FOR NON-LABORATORY SETTINGS: In a non-laboratory setting such as an endoscopy suite, at least one per lot of *hpf*ast, a quality control test should be performed. A single drop of the urease-glycerol preparation should be inserted in an appropriate number of *hpf*ast samples. A single drop should turn a test positive within one hour.

8.6.3 FAIL SAFE TESTING FOR CLIA WAIVER: In a non-laboratory setting the positive control solution should also be used to verify the test integrity of any negative test. Place a single drop of the urease in the *hpf*ast and re-seal the test. The urease should turn the test positive within one hour.
⚠ CAUTION: *Immediately mark the test so it cannot be confused as a true positive patient specimen! Discard all positive controls immediately after reading the test.*

8.6.4 POSITIVE CONTROLS FOR LABORATORIES: In a laboratory setting, a viable culture of *Helicobacter pylori* or another strong urease producing organism such as *Proteus mirabilis* available from the American Type Culture Collection (tel. # 301-881-2600, ATCC strain # 33583) may be substituted for the commercially prepared urease. Since *hpf*ast is highly specific for *H. pylori*, a large amount of the substituted organism (at least one million organisms) may be necessary to turn the test positive. The *hpf*ast should turn positive in approximately 1 hour.

NOTE: Due to the high specificity of hpfast and the difficulties associated with maintaining and quantifying viable cultures, use of the provided urease-glycerol preparation is highly recommended as the positive control.

8.6.5 NEGATIVE CONTROLS: The color of a negative test (yellow/peach color) is indicated on the *hpf*ast label. A non-urease producing organism or a single drop of sterile saline solution may be used as a negative control.

Should any of the controls react incorrectly, contact GI Supply at 1-800-451-5797.

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9. LIMITACIONES AND WARNINGS

9.1 Falso Positivo

Studies indicate that the urease of other organisms such as *Proteus* or *Pseudomonas* is not as active in the acidic pH range where the *hpf*ast reaction begins. This feature tends to separate out *H. pylori* urease activity from that in other bacteria. Patients with complete achlorhydria (pernicious anemia, previous gastric surgery and long term acid suppression by drugs) could have gastric bacterial growth of other urease containing bacteria. In general, these bacteria have less urease and do not produce a rapid or deep blue positive test.

9.2 Falso Negativo

Several factors may produce a false negative result.

a. Patchy Distribution of Infection - *H. pylori* has been shown to occur in a patchy distribution so that two antral biopsies are recommended, usually in the immediate prepyloric area and on the lesser curvature of the proximal antrum.

b. Intestinal Metaplasia - *H. pylori* does not colonize intestinal metaplasia so that if this mucosal change is extensive, a false result may occur. Again, multiple biopsies are recommended.

c. Antibiotics - *H. pylori* is sensitive to many antimicrobials, including bismuth. If patients have been on antimicrobials in the recent past, the organism may be suppressed though usually not eliminated. In this situation, an inadequate number of urease containing bacteria may be obtained. If recent antibiotics have been used, multiple biopsies should be taken.

d. Omeprazole - This proton pump inhibitor has been shown to have rapid bacteriocidal effects against *H. pylori*. Ideally, the drug should be discontinued for two weeks before biopsy.

9.3 In difficult or questionable cases of infection, additional testing by histology, serology, or labeled carbon breath tests (when available) should be considered. *hpf*ast is intended for in-vitro diagnostic use only.

10. EXPECTED VALUES

For the symptomatic population in the *hpf*ast clinical study, 111 out of 470 patients or 23.6% were positive by the *hpf*ast test criteria.

11. PERFORMANCE CHARACTERISTICS

In a study of 470 patients with possible peptic ulcer disease, 122 demonstrated positive histological criteria (acute gastritis and presence of *H. pylori* organisms upon staining) for *H. pylori* gastritis and 348 did not demonstrate the criteria. *hpf*ast agreed with the positive histological diagnosis in 111 cases and agreed with the negative diagnosis in 342 of the cases. *This data provides a relative sensitivity of 90.97% and a relative specificity of 98.27%.*

| | | |
|--------------------------|----------------------|----------------------|
| <i>hpf</i> ast Positives | True + 111 | False + 6 |
| <i>hpf</i> ast Negatives | False — 11 | True — 342 |

Histology Positives = 122


Histology Negatives = 348

*hpf*ast Sensitivity = 90.97%

*hpf*ast Specificity = 98.27%

*hpf*ast diagnosed 76.6% of the positives by 30 minutes, 89.2% were positive by 60 minutes, and 91.9% by 4 hours. There were no false positives when diagnosed within the first hour of the test, only one false positive when diagnosed within four hours, and six false positives when diagnosed at 24 hours.

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