5 1/2" (+/- 1/16")

ENDOSCOPIC TATTOO INSTRUCTIONS FOR USE

DESCRIPTION

SPOT® Endoscopic Tattoo (SPOT) is a specially formulated, biocompatible suspension designed for endoscopically tattooing or marking lesions in the gastrointestinal (GI) tract. The suspension contains highly purified, very fine carbon particles which remain in the submucosal injection site creating a permanent tattoo. Each steam sterilized, non-pyrogenic syringe contains 8ml of the following ingredients:

- High Purity Carbon Black
- Polysorbate 80
- Glycerol
- Water
- Simethicone
- Benzyl Alcohol

INDICATION FOR USE

SPOT is indicated for endoscopically marking lesions in the GI tract when the endoscopist anticipates the lesion will require surgical removal within 30 days.

CONTRAINDICATIONS

SPOT is contraindicated in patients with a known hypersensitivity to any of the ingredients in the product. **Not For Pediatric Use.**

WARNING AND PRECAUTIONS

Some settling of the carbon pigment is normal. The syringe should be shaken vigorously for 15-20 seconds prior to loading the injection needle.

SPOT is a sterile single-patient-use product. Sterile unless damaged. Leakage of product may be evidence of damage. Do not use if damaged. Discard any unused product after the syringe has been opened. To avoid cross contamination, each syringe may only be used on one (1) patient and cannot be reused. The suspension may stain clothing or porous materials.

When using an injection needle in the GI tract, the needle should be inserted into the GI wall at an angle so that the beveled tip of the needle is entirely beneath the mucosa. If the needle is inserted perpendicular to the mucosa, there is a risk of penetrating the serosa and injecting the tattoo suspension directly into the peritoneal cavity.

SIDE EFFECTS

Any inadvertent injection of tattoo suspension into the intraperitoneal space should be monitored closely for adverse reactions.

DIRECTIONS FOR USE

Open the tray and remove the syringe.

IMPORTANT: The SPOT syringe should be shaken vigorously for 15-20 seconds prior to loading the injection needle. Examine the syringe to verify that the pigment is fully suspended. Then, using aseptic technique, attach the syringe to the Luer fitting on the needle. A 23 or 25 gauge injection needle with a needle length of 4mm or less is recommended for this procedure. **Prime needle with SPOT prior to injection.** With the needle retracted, insert the catheter through the biopsy working channel of the endoscope.

When the needle is properly positioned, insert the tip 30°-45° tangentially into the submucosal space of the GI wall. **CAUTION: Do not insert the needle perpendicular to the GI wall as this may lead to perforation of the serosa and injection of the tattoo suspension directly into the peritoneal cavity!**

Inject 0.5-0.75 ml of the tattoo into each injection site. This volume typically produces a mark on the GI wall measuring 1-2 cm in diameter. When properly placed, there is an immediate black blush of the tattoo at the injection site. Observe for leakage of the tattoo into the lumen. If this occurs, stop the injection and cautiously insert the needle further into the submucosa. Consider using the saline test bleb method which first injects 1 ml of saline to create a mound or bleb prior to injecting SPOT. The saline test bleb method confirms placement in the submucosal layer.

Three to four tattoo injections are recommended to create a circumferential pattern around the lesion. This technique assures that at least one of the marks can be observed during surgery from the serosal side of the GI tract.

The maximum recommended dose is 8 ml per patient.

STORAGE AND SHELF LIFE

Store the product at room temperature. The product has a shelf life of two years. The expiration date and the lot number are printed on the syringe label. Store syringe with the syringe tip pointed upward. **Do not store with the syringe tip pointing down.** Storing the syringe tip down may allow a portion of the carbon pigment to settle into the narrow syringe outlet and prevent it from being re-suspended.

CUSTOMER SERVICE

GI Supply, 5069 Ritter Road, Suite 104, Mechanicsburg, PA 17055 USA • (800) 451-5797 • 8:30am and 5:00pm EST
Monday-Friday www.gi-supply.com • orders@gi-supply.com

Sterile Fluid Path

Use by (Expiration Date)

Lot number

Consult Instructions for Use

Caution

Non Pyrogenic

Manufactured By

Caution: Federal law restricts this device to sale by or on the order of a physician

Lower Limit Temperature

Keep Dry

Do Not Use if Package is Damaged

SPOT is a registered trademark of GI Supply, Inc.
©2020 GI Supply, Inc. All rights reserved.

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©2020 GI Supply, Inc I Todos los